


































**3**

### Assay Data Sheet - Expected Ranges BioProfile® FLEX Controls - LEVEL 3

**3**

Lot Number:  427317				Exp. Date: 2016-03	
Expected Range				+/- 3 SD Expected Range	
Parameter	Units	Lower Limit	Upper Limit	Lower Limit	Upper Limit
pH	-	 7.326	 7.426	 7.301	 7.451
PO <sub>2</sub>	mm Hg	 139.9	 168.5	 132.8	 175.6
PCO <sub>2</sub>	mm Hg	 14.9	 21.3	 13.3	 22.9
Gluc	g/L	 12.75	 17.25	 11.62	 18.38
Lac	g/L	 4.25	 5.75	 3.87	 6.13
NH <sub>4</sub> <sup>+</sup>	mmol/L	 8.24	 11.24	 7.49	 11.99
Na <sup>+</sup>	mmol/L	 166.7	 203.3		
K <sup>+</sup>	mmol/L	 5.55	 6.85		
Ca <sup>++</sup>	mmol/L	 2.25	 2.85		
Osmolality	mOsm/kg H <sub>2</sub> O	 601	 617		

*nova*  
biomedical

200 Prospect Street • Waltham, MA U.S.A. 02454-9141 • (781) 864-0800

LFN 08241A 2007-08